



Limited liability company

«Biolabmix»

TIN 5408278957 CAT 540801001

630090, Novosibirsk obl., Novosibirsk,

st. Injenernaya, building № 28

Tel/Fax: +7(383)363-51-91, Tel: +7(383)363-22-40

E-mail: sales@biolabmix.ru

BioMaster HS-*Taq* PCR Kit

Cat. number KH017-500, KH017-2250

Product description:

HS-*Taq* PCR Kit contains recombinant HS-*Taq* DNA polymerase and the solutions containing all the essential components to perform conventional hot-start PCR, excluding DNA template and primers. The kit includes: HS-*Taq* DNA polymerase solution (5 U/μl), 5× PCR buffer, 50 mM MgCl₂, 50× dNTP mix and 6× loading buffer. All of the included reagents are high quality and optimized for PCR.

HS-*Taq* DNA polymerase is a recombinant *Taq* DNA polymerase inactivated by specific monoclonal antibodies. HS-*Taq* DNA polymerase is inactive at temperatures lower than 70 °C. It allows to avoid the formation of non-specific products and primer-dimers at low temperatures during sample preparation. Activation of HS-*Taq* DNA polymerase occurs during the first cycle of PCR amplification after short 5 min incubation at 95 °C. The recombinant form of the enzyme has 5'-3' DNA-dependent DNA polymerase activity and 5'-3' exonuclease activity of native *Taq* DNA polymerase from *Thermus aquaticus*. The extension rate of *Taq* DNA polymerase depends on the complexity of DNA template and is approximately 2 kbp/min. Recombinant *Taq* DNA polymerase is ideal for conventional PCR of templates up to 5 kbp.

5× PCR buffer (**does not contain MgCl₂**) is optimized for efficient and reproducible PCR. The buffer components increase half-life and processivity of HS-*Taq* DNA polymerase by enhancing its stability during PCR. The buffer is chemically stable, inert and its composition does not affect optimal annealing temperature or the parameters of template melting. Additional solutions of 50 mM MgCl₂ and 50× dNTP mix included in the kit allow optimizing PCR conditions for individual primer pairs, while 6× loading buffer eases sample preparation for analysis and allows monitoring the progress of electrophoresis.

Product composition

| Cat. # | HS- <i>Taq</i> DNA polymerase, 5 U/μl* | 5× PCR buffer | 50 mM MgCl ₂ | 50× dNTP mix (10 mM each) | 6× loading buffer | Activity units (U) |
|------------|--|---------------|-------------------------|---------------------------|-------------------|--------------------|
| KH017-500 | 1 × 100 μl | 2 × 1.5 ml | 1 × 1 ml | 2 × 200 μl | 1 × 1.75 ml | 500 |
| KH017-2250 | 3 × 150 μl | 8 × 1.5 ml | 2 × 1 ml | 4 × 400 μl | 3 × 1.75 ml | 2250 |

* One activity unit is the amount of the enzyme required to catalyze incorporation of 10 nmol of dNTP into an acid-insoluble product in 30 min at 74 °C. Reaction conditions: 50 mM Tris-HCl (pH 9.0 at 25 °C), 50 mM NaCl, 10 mM MgCl₂, 200 mM dATP, 200 mM dCTP, 200 mM dGTP, 50 mM [³H] dTTP, 0.25 mg/ml of activated calf thymus DNA.

Storage buffer:

50 mM Tris-HCl (pH 8.0 at 25 °C), 50 mM NaCl, 0.1 mM EDTA, 1 mM dithiothreitol, 50% (v/v) glycerin, and 1% (v/v) Triton X-100.

5× PCR buffer:

50 mM Tris-HCl (pH 8.5 at 25 °C), 250 mM KCl, 0.5% (v/v) Tween 20, stabilizers of *Taq* DNA polymerase.

Applications:

- Hot-start PCR;
- High-throughput PCR;
- Conventional PCR with high reproducibility;
- Synthesis of PCR products for TA cloning;
- RT-PCR.

Limits of use

- Not recommended to use for amplicons of >5 kbp.

Inhibition and inactivation

Inhibitors: ionic detergents (sodium deoxycholate, sarkosyl and sodium dodecyl sulfate (SDS) at concentrations above 0.06, 0.02 and 0.01%, respectively).

Inactivated by phenol/chloroform extraction.

Amplification protocol

Prepare several parallel reactions to minimize possible pipetting error: mix water, buffer, dNTP mix, primers and HS- *Taq* DNA polymerase. Prepare reaction solution by multiplying the total volume by the number of reactions required plus one additional reaction to cover pipetting losses. Aliquot PCR solution to individual PRC tubes, then add DNA template.

1. Thaw the reaction mix and mix carefully.

Note: in case of precipitate formation, heat the tube at 50 °C and stir till complete dissolution.

2. Add the following components into thin-wall PCR tubes considering that the final volume of a reaction mixture is 50 µl:

| Component | Volume | Final concentration |
|---------------------------------------|-------------|---------------------|
| 5× PCR buffer | 10 µl | 1× |
| 50× dNTP mix | 1 µl | 0.2 mM of each |
| 50 mM MgCl ₂ * | variable | 1-5 mM |
| Forward primer | variable | 0,1 – 300 nM |
| Reverse primer | variable | 0,1 – 300 nM |
| DNA template | variable | 10 pg – 1 µg |
| HS- <i>Taq</i> DNA polymerase, 5 U/µl | variable | 1-5 U |
| Sterile water | up to 50 µl | |

3. Carefully vortex and collect all droplets from the tube walls by brief centrifuging.

Note: in case if a thermal cycler is not equipped with a heated lid, add a droplet (25-35 µl) of mineral oil in each tube.

4. Perform PCR using conditions recommended below:

| Step | Temperature, °C | Incubation time | Number of cycles |
|--------------------------|-----------------------------|-----------------|------------------|
| Preliminary denaturation | 95 | 5 min | 1 |
| Denaturation | 95 | 5 – 10 sec | |
| Annealing | 50 – 68 (T _m -5) | 5 - 15 sec | 25 - 45 |
| Elongation | 72 | 1-2 min/kbp | |
| Final elongation | 72 | 5 – 15 min | 1 |

T_m – template/primer melting temperature, depends on primer structure. The following formula can be used for approximate estimation of T_m: $T_m (°C) = 2 \times (A+T) + 4 \times (G+C)$.

5. After performing PCR, analyze amplification products by gel electrophoresis. Samples are loaded on gel without additional loading buffer.

Note: we recommend using 1x TAE buffer with ethidium bromide for visualizing PCR products by gel electrophoresis.

Note: mobility of dyes in 0.5 – 1.5% agarose gel:

| xylene cyanol | bromphenol blue | Orange G | tartrazine |
|-----------------|-----------------|----------|------------|
| 10000 – 4000 bp | 500-400 bp | <100 bp | <20 bp |

Storage conditions: in a place protected from light at +25 ° C - 7 days; at +4 ° C - 4 months; at -20 ° C - 12 months; not more than 50 thawing-freezing cycles.

Transportation: Transport in thermocontainers with cooling elements; the ambient temperature increment to the room temperature during the transportation up to 10 days is allowed.