



Limited liability company  
**«Biolabmix»**  
TIN 5408278957 CAT 540801001  
630090, Novosibirsk obl., Novosibirsk,  
st. Injenernaya, building № 28  
Tel/Fax: +7(383)363-51-91, Tel: +7(383)363-22-40  
E-mail: sales@biolabmix.ru

## BioMaster RT-PCR – Standard (2×)

Cat. Number: RM02-40, RM02-200

### Product description:

**BioMaster RT-PCR – Standard (2×)** was designed for one-tube reverse transcription and polymerase chain reaction (RT-PCR). The kit includes **2× buffer for RT-PCR-Standard** containing all the necessary components (excluding DNA template and primers); enzyme mix **BioMaster-mix** and **DEPC-treated water, DMSO** and **6×loading buffer**.

**BioMaster-mix** contains the optimal ratio of M-MuLV –RH to *HS-Taq* DNA polymerase for both reactions.

M-MuLV –RH is a genetically modified Moloney Murine Leukemia Virus reverse transcriptase (M-MuLV). The enzyme exerts RNA- and DNA-dependent polymerase activity but lacks RNase H activity. M-MuLV –RH reverse transcriptase exhibits improved thermal stability and remains active up to 50°C.

*HS-Taq* DNA polymerase is a recombinant *Taq* DNA polymerase inactivated by specific monoclonal antibodies. The enzyme is inactive at temperatures up to 70 °C, it is activated at the first PCR cycle during a short 5-min incubation at 95 °C. Recombinant *HS-Taq* DNA polymerase catalyzes 5'→3' synthesis of DNA and possesses 5'→3' exonuclease activity of the native *Taq* DNA polymerase from *Thermus aquaticus*. Recombinant *HS-Taq* DNA polymerase is ideal for conventional PCR of the templates of up to 5 kbp.

The buffer is optimized for both efficient RT and PCR.

### Product composition:

Cat. #	2× buffer for RT-PCR-Standard	BioMaster-mix	DMSO	DEPC-treated water	6×loading buffer
RM02-40	2 × 0.5 ml	1 × 80 µl	1 × 0.1 ml	2 × 0.5 ml	1 × 0.1 ml
RM02-200	4 × 1.25 ml	1 × 400 µl	1 × 0.2 ml	3 × 1.8 ml	1 × 0.5 ml

### 2× buffer for RT-PCR-Standard contains:

100 mM Tris-HCl (pH 8.3 at 25 °C), 150 mM KCl, 0.6 mM each deoxynucleoside triphosphate, 6 mM MgCl<sub>2</sub>, 8 mM DTT, enzyme stabilizers and enhancers.

### BioMaster-mix contains:

50 mM Tris-HCl (pH 8.0 at 25 °C), 100 mM NaCl, 1 mM EDTA, 5 mM DTT, 50 % (v/v) glycerol and 0.1 % (v/v) NP-40, M-MuLV –RH reverse transcriptase and *HS-Taq* DNA polymerase.

### Applications:

- Gene expression analysis;
- One-step conventional RT-PCR.

## Reaction mix features

- The reaction mix is optimized for the specific and effective performance of M-MuLV –RH reverse transcriptase and *HS-Taq* DNA polymerase;
- Allows long-term storage (storage of **BioMaster RT-PCR – Standard (2×)** at room temperature for 2 days and/or multiple thawing-freezing cycles do not affect RT-PCR efficacy);

## Benefits of use

- High specificity;
- High sensitivity;
- Convenient and easy-to-use;
- Low pipetting error and low risk of cross-contamination;
- Standardized conditions of the same-type reactions (reduced pipetting error during mixing PCR components in a series of experiments);
- PCR products can be further subjected to TA cloning due to deoxyadenosine overhangs at the ends of amplified DNA fragments;
- Shorter step of sample preparation for the analysis of PCR results.

## Limits of use

- Not recommended for amplicons of > 5 kbp.

## Protocol

Before starting to work, we recommend to get acquainted with the protocol and recommendations presented at our site: <http://biolabmix.ru/catalog>

1. Thaw **2× buffer for RT-PCR-Standard** and vortex thoroughly.
2. Place the thin-wall tubes in ice and add the following components considering the final volume of a reaction mixture equal to 50 µl:

Component	Volume	Final concentration
2× mix for RT-PCR-Standard	25 µl	1×
BioMaster-mix	2 µl	
Forward primer	variable	0.1 – 500 nM
Reverse primer	variable	0.1 – 500 nM
RNA template	variable	1 pg – 1 µg
Sterile water	up to 50 µl	

**Note:** in case of amplification of templates with complicated spatial structure, DMSO can be added in the amount of 1 to 5% of the final volume of the reaction solution. Change in  $T_m$  of the primers should be taken into account when selecting the amplification program.

**Note:** the volume of **BioMaster-mix** can be varied in the range of 1 to 3 µl per 50 µl reaction depending on the gene copy number and complexity.

3. Carefully vortex and discard the droplets by centrifugation.

**Note:** in case of using the thermal cycler without a heating lid, add a drop of mineral oil (25-35  $\mu$ l) to each tube.

4. Perform PCR using recommended conditions:

Step	Temperature, °C	Incubation time	Number of cycles
Reverse transcription	45	30 min	1
Preliminary denaturation	95	5 min	1
Denaturation	95	5 – 15 sec	25 – 45
Annealing	50 – 68 ( $T_m-5$ )	5 – 20 sec	
Elongation	72	0.5-1 min/kbp	
Final elongation	72	5 – 15 min	1

$T_m$ : template-primer duplex melting temperature, it depends on the primer structure.

The following formula can be used for  $T_m$  estimation:

$$T_m (^{\circ}\text{C}) = 2 \times (\text{A}+\text{T}) + 4 \times (\text{G}+\text{C}).$$

5. After conducting PCR, analyze amplification products by gel electrophoresis. No loading buffer is required.

**Note:** we recommend using 1xTAE buffer with ethidium bromide for separation of amplification products by gel electrophoresis.

**Note:** mobility of dyes in 0.5 – 1.5% agarose gel:

xylene cyanol	bromphenol blue	Orange G	tartrazine
10000 – 4000 bp	500-400 bp	<100 bp	<20 bp

### Optimization of reaction conditions

1. The reaction volume can be varied in the range of 10 to 50  $\mu$ l with proportional change in the amount of all components.
2. When using a template containing GC-rich regions and regions with complicated spatial structure, the temperature can be increased to 50  $^{\circ}$ C, and/or reagents facilitating melting of the secondary structure of the nucleic acids (e.g. DMSO) can be added.

**Storage conditions:** in a place protected from light at +4  $^{\circ}$ C – 1 month; at -20  $^{\circ}$ C - 1 year; not more than 30 thawing-freezing cycles.

**Transportation** Transport in thermocontainers with cooling elements; the ambient temperature increment to the room temperature during the transportation up to 10 days is allowed.