



Limited liability company  
«Biolabmix»  
TIN 5408278957 CAT 540801001  
630090, Novosibirsk obl., Novosibirsk,  
st. Injenernaya, building № 28  
Tel/Fax: +7(383)363-51-91, Tel: +7(383)363-22-40  
E-mail: sales@biolabmix.ru

## Thermolabile alkaline phosphatase

Cat. Number: E-12005, E-12050

### Enzyme description

This product is a recombinant enzyme – alkaline phosphatase of the gram-negative bacterium *Vibrio splendidus*. The enzyme has a molecular weight of ~58 kDa and exhibits catalytic activity in its homodimeric form [1]. The temperature range of the enzyme is from 15 to 37°C, but for prolonged reactions it is recommended to use 15°C as the most optimal for enzyme stability [2]. Thermolabile alkaline phosphatase removes phosphates from the 5' and 3' ends of DNA and RNA and can be used instead of Antarctic phosphatase.

### Application

Cloning of restriction fragments.  
mRNA synthesis.

### Source

Thermolabile alkaline phosphatase is isolated from an *E. coli* strain containing a plasmid with a cloned *Vibrio splendidus* enzyme gene.

### Activity Units

One unit of activity hydrolyzes 1 µmol of para-nitrophenylphosphate (pNPP) in 0.5 ml of the reaction mixture in 15 min at 25°C, in a standard reaction buffer (10 mM Tris-HCl (pH 9.0 at 25°C); 10 mM MgCl<sub>2</sub>; 100 mM NaCl; 1 mM DTT, 10 mM pNPP).

**Enzyme concentration and packaging:** 5000 U/ml.

Cat. No.	Product Name	Quantity	Volume
E-12005	Thermolabile alkaline	500 U	100 µl
E-12050	phosphatase	5000 U	1000 µl

### Storage Buffer

25 mM Tris (pH 8.0 at 25°C), 2 mM MgCl<sub>2</sub>, 0.02 mM ZnCl<sub>2</sub>, 50% glycerol

### Quality Control

Each batch of enzyme is tested for enzyme activity, electrophoretic purity in SDS-PAGE, and nonspecific proteolytic activity.

### **10x Standard reaction buffer**

**100 mM Tris-HCl (pH 9.0 at 25°C), 100 mM MgCl<sub>2</sub>, 1 M NaCl, 10 mM DTT (1 ml 10x buffer supplied with enzyme).**

### **Typical protocol for dephosphorylation of hydrolyzed plasmids**

1. Mix the following components in a test tube:
  - 2 µl of 10x standard reaction buffer;
  - 1 µg of hydrolyzed plasmid DNA (4–8 kb);
  - up to 19 µl of nuclease-free water;
  - 1 µl (5 U) thermolabile alkaline phosphatase.
2. Incubate the reaction mixture at 37°C for one hour.
3. To inactivate the enzyme, heat the reaction mixture at 65°C for 20 minutes.
4. The processed plasmid DNA can be used in further work. If necessary, DNA can be further purified using well-known methods, for example, a Kit for isolating DNA and RNA from reaction mixtures (cat. no. DR-10, DR-50, DR-250).

### **Enzyme inactivation**

Incubation for 20 minutes at 65°C.

### **Storage and transportation conditions**

Store at -20°C.

Transportation is allowed at a temperature not exceeding +8°C for two days.

### **References**

1. Helland, R., Larsen, R. L., & Ásgeirsson, B. (2009). The 1.4 Å crystal structure of the large and cold-active *Vibrio* sp. alkaline phosphatase. *Biochimica et Biophysica Acta (BBA)-Proteins and Proteomics*, 1794(2), 297–308.
2. Hauksson, J. B., Andrésón, Ó. S., & Ásgeirsson, B. (2000). Heat-labile bacterial alkaline phosphatase from a marine *Vibrio* sp. *Enzyme and microbial technology*, 27(1–2), 66–73.