



Limited liability company

«Biolabmix»

TIN 5408278957 CAT 540801001

630090, Novosibirsk obl., Novosibirsk,

st. Injenernaya, building № 28

Tel/Fax: +7(383)363-51-91, Tel: +7(383)363-22-40

E-mail: sales@biolabmix.ru

Bst DNA-polymerase

Cat. Number: E-10002, E-10010

Product description

DNA-polymerase LF Bst is a large Bst polymerase fragment of *Bacillus stearothermophilus*. Enzyme contains histidine mark on its C-terminus and has molecular mass of 68,9 kDa. The enzyme is high-processive and catalyses 5'→3' DNA synthesis. Fragment lacks 5'→3' and 3'→5'-exonuclease activity as well as 5'-3' displacement activity.

Optimal activity conditions for the enzyme are 65 °C and pH 8,9.

Application:

- loop isothermal amplification (LAMP)
- Whole-genome sequencing

Source:

Bst DNA-polymerase was obtained from *E.coli* strain, carrying plasmid with cloned full-size gene of *Bacillus stearothermophilus* DNA-polymerase I large fragment.

Activity units:

One activity unit corresponds to the enzyme amount required for inclusion of 10 nmoles dNTPs in non-acid-soluble DNA fraction in 0 min at 65°C.

Quality control:

Every enzyme batch is tested for endonuclease and non-specific exonuclease activity absence, enzyme sensitivity.

Enzyme concentration: 10 UA/μl.

Storage buffer: to the amount of enzyme 20 mM KCl, 1 mM DTT, 0.1 mM EDTA, 0.1% Triton X-100, 50% glycerin.

10x LAMP guffer (Cat №: SPO30, supplied separately):

300 mM Tris-HCl (pH 8.9), 50 mM (NH₄)₂SO₄, 0.5 mg/ml BSA, 2.0% Tween 20.

Typical LAMP assay:

Reaction mixture contains:

Component	Final concentration	Volume, μ l
10 \times LAMP buffer (Cat N $^{\circ}$: SP030)	1 \times	5
100 mM MgSO ₄	4–10 mM (6 mM)*	3*
10 mM dNTP Mix (Cat N $^{\circ}$: NM10)	1.4 mM each	7
16 μ M FIP/BIP Primers	1.6 μ M	5
2 μ M F3/B3 Primers	0.2 μ M	5
8 μ M LoopF/B Primers	0.8 μ M	5
10 U/ μ l <i>Bst</i> DNA-polymerase	0.008 – 0.2 10 U/ μ l (0.033 U/ μ l)*	0.2*
DNA-matrix	> 10 copies per reaction	variable
Sterile water (Cat. N $^{\circ}$: SP010)		up to 50
Total reaction volume	50 μ l	

* recommended concentration

To perform amplification specialised isothermal reaction amplifiers as well as PCR amplifiers. It is recommended to perform amplification reaction at the 60 – 65 °C temperature range with the duration of 20–30 minutes. For the correct results interpretation it is recommended to always set a control reaction without the DNA-matrix.

Transportation: in thermostated containers with cooling elements, tolerating temperature increment up to environment temperature, if transported in 10 days.

Storage conditions and period: 1 year at -20 °C.