




Limited liability company
«Biolabmix»
TIN 5408278957 CAT 540801001
630090, Novosibirsk obl., Novosibirsk,
st. Injenernaya, building № 28
Tel/Fax: +7(383)363-51-91, Tel: +7(383)363-22-40
E-mail: sales@biolabmix.ru

N6-methyladenosine-5'-triphosphate (N6-Me-ATP)

Catalog number: TNA-0050, TNA-0500



N6-Methyl-ATP, m6ATP

Molecular formula (acid): C₁₁H₁₈N₅O₁₃P₃

Molecular weight (acid): 521,21 g/mol

Description:

The product is a sterile solution of 100 mM of N6-Me-ATP ammonium salt in water. The nucleotide has not less than 96% purity confirmed by HPLC, is tested for the presence of endo- and exonuclease activity and free of trace contaminating DNase and RNase, functionally tested *in vitro* in the transcription reaction.

FOR RESEARCH USE ONLY.

Packing size

Cat. number	Product name	Composition	Quantity
TNA-0050	N6-Me-ATP	100 mM N6-Me-ATP ammonium salt in water	50 µl
TNA-0500	N6-Me-ATP	100 mM N6-Me-ATP ammonium salt in water	500 µl

Storage and transportation conditions:

Long term storage at -20°C. It is transported in a frozen form, the solution may be defrosted and transported at ambient temperature.

Shelf life: 12 months.